


1



2



3

$$\Delta = (Pwr * dec)/10$$

Where:

- Δ = prism diopters
- Pwr = Dioptic value
- dec = decentration (in mm)

4

Review of ANSI:

- What does the ANSI acronym stand for?
- What does purpose does ANSI serve?
- IN Opticianry, What are the ANSI standards we use most often?
-

5

Review of ANSI:

ANSI is a VOLUNTARY industry standard.

Discuss.

6

Review of ANSI:

ANSI Z80.1 – 2020 Prescription Ophthalmic Lenses, Recommendations

State Opticians laws in licensed states (CT and MA specifically)

ANSI Z87.1 – 2020 Occupational and Educational Personal Eye and Face Protection Devices
OSHA Federal Agency

7

HIGHLY RECOMMEND YOU PURCHASE ENTIRE STANDARDS

[https://webstore.ansi.org/Standards/VC%20\(ASC%20Z80\)/ANSIZ802020-2436966](https://webstore.ansi.org/Standards/VC%20(ASC%20Z80)/ANSIZ802020-2436966)

8

Most folks keep a "cheat sheet" for the ANSI Z80.1 next to their lensometer for reference...


Make sure it is the latest, 2020

ALSO make sure you know how to read it.

https://thevisioncouncil.org/sites/default/files/assets/media/ANSI_Z80-1-2020_QuickReferenceGuide.pdf


due to potential copyright we will not distribute this

9





13

13




14

14




15

15



16



17




Just to confuse the matter, there are even MORE!!!!


Today's lab has a few different reticle types

18







19



20



21



22

Prism notation

Horiz amount and direction with Vertical Amount and direction

Or


Base at a specified Angle

23


23

Directions of Prism in a Lensmeter


Reticle viewed as though looking through eyepiece



24




25



26



27



28

Given the following Rx:

OD -1.00 sph prism 3 Δ @135
 OS -2.00 sph prism 3 Δ @ 315
 Add +2.25 OU

The prism in the OD would best be described as:

- Up and Out
- Up and In
- Down and Out
- Down and IN

The Prism in the OS would best be described as:

- Up and Out
- Up and In
- Down and Out
- Down and IN

29

29

ANSI standards and applications:

- Check Dist Power tolerance
- Check Cyl Power Tolerance
- Check Axis Tolerance
- Check ADD

- Determine vertical imbalance
- Determine Horizontal imbalance

- Seg vert and horizontal or Fitting Cross Vert and Horiz


30

30


ANSI standards and applications:

- Also don't forget
 - BC (really?)
 - CT
 - Segment Size and tilt

31

<p>ORDERED</p> <p>-2.00 -0.25 x 180</p> <ul style="list-style-type: none"> • Step 1 check SPH power <ul style="list-style-type: none"> ◦ Accept ANYTHING within 0.13 ◦ OK from -1.87 to -2.13 • Step 2 check CYL power <ul style="list-style-type: none"> ◦ Accept anything within 0.13 ◦ OK from -0.12 to -0.38 • Step 3 check Cyl AXIS <ul style="list-style-type: none"> ◦ Accept anything within 14 degrees ◦ 166 to 014 	<p>RECEIVED</p> <p>-2.00 -0.25 x 005</p> 
--	---

32

<p>ORDERED</p> <p>-2.00 -0.25 x 180</p> <ul style="list-style-type: none"> • Step 1 check SPH power <ul style="list-style-type: none"> ◦ Accept ANYTHING within 0.13 ◦ OK from -1.87 to -2.13 • Step 2 check CYL power <ul style="list-style-type: none"> ◦ Accept anything within 0.13 ◦ OK from -0.12 to -0.38 • Step 3 check Cyl AXIS <ul style="list-style-type: none"> ◦ Accept anything within 14 degrees ◦ 166 to 014 	<p>RECEIVED</p> <p>-1.88 -0.33 x 168</p> 
--	--

33

<p>ORDERED</p> <p>-2.00 -0.25 x 180</p> <ul style="list-style-type: none"> · Step 1 check SPH power <ul style="list-style-type: none"> ○ Accept ANYTHING within 0.13 ○ OK from -1.87 to -2.13 · Step 2 check CYL power <ul style="list-style-type: none"> ○ Accept anything withing 0.13 ○ OK from -0.12 to -0.38 · Step 3 check Cyl AXIS <ul style="list-style-type: none"> ○ Accept anything within 14 degrees ○ 166 to 014 	<p>RECEIVED</p> <p>-2.12 -0.15 x 013</p> <div style="background-color: #ADD8E6; width: 150px; height: 100px; margin: 10px 0;"></div> <p>ansi, so</p>
---	--

34

<p>ORDERED</p> <p>-2.00 -0.25 x 180</p> <ul style="list-style-type: none"> · Step 1 check SPH power <ul style="list-style-type: none"> ○ Accept ANYTHING within 0.13 ○ OK from -1.87 to -2.13 · Step 2 check CYL power <ul style="list-style-type: none"> ○ Accept anything withing 0.13 ○ OK from -0.12 to -0.38 · Step 3 check Cyl AXIS <ul style="list-style-type: none"> ○ Accept anything within 14 degrees ○ 166 to 014 	<p>RECEIVED</p> <p>-2.00 -0.25 x 015</p> <div style="background-color: #ADD8E6; width: 150px; height: 100px; margin: 10px 0;"></div>
---	--

35

<p>ORDERED</p> <p>-10.00 -3.00 x 0180</p> <ol style="list-style-type: none"> 1. check SPH power <ul style="list-style-type: none"> ○ Accept ANYTHING within 2% ○ Accept anything withing 0.2 ○ Anything between -9.80 and -10.20 is within tolerance 2. check CYL power <ul style="list-style-type: none"> ○ Accept anything within 0.15 ○ Accept anything between -2.85 and -3.15 3. check Cyl AXIS <ul style="list-style-type: none"> ○ Accept anything within 2 degrees ○ 178-002 	<p>If RECEIVED:</p> <p>-9.85 -3.10 x 178</p> <div style="background-color: #ADD8E6; width: 150px; height: 20px; margin: 5px 0;"></div> <p>If Received</p> <p>-10.00 -3.00 x 003</p> <div style="background-color: #ADD8E6; width: 150px; height: 20px; margin: 5px 0;"></div>
---	---

36

ORDERED


2% = power * 2, and move decimal to the left 2 spots

- 2% of 10 = 0.2
- 2% of 9 = .18
- 2% of 14 = .28

37

You order the following pair of eyewear:

38



39
